

English Version

Plant biostimulants - Determination of inorganic arsenic

Biostimulants des végétaux - Dosage de l'arsenic
inorganique

Pflanzen-Biostimulanzien - Bestimmung von
anorganischem Arsen

This draft European Standard is submitted to CEN members for enquiry. It has been drawn up by the Technical Committee CEN/TC 455.

If this draft becomes a European Standard, CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

This draft European Standard was established by CEN in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Türkiye and United Kingdom.

Recipients of this draft are invited to submit, with their comments, notification of any relevant patent rights of which they are aware and to provide supporting documentation.

Warning : This document is not a European Standard. It is distributed for review and comments. It is subject to change without notice and shall not be referred to as a European Standard.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

| Contents | Page |
|--|-------------|
| European foreword | 3 |
| Introduction | 4 |
| 1 Scope..... | 5 |
| 2 Normative references..... | 5 |
| 3 Terms and definitions | 5 |
| 4 Principle | 5 |
| 5 Reagents | 6 |
| 6 Apparatus | 8 |
| 7 Sampling..... | 8 |
| 8 Procedure | 9 |
| 8.1 Sample preparation | 9 |
| 8.2 Water bath extraction | 9 |
| 8.3 Preparation of the calibration solutions | 9 |
| 8.4 Measurement..... | 9 |
| 8.4.1 General instrumental conditions | 9 |
| 8.4.2 HPLC-ICP/MS measurement of inorganic arsenic..... | 10 |
| 9 Calculation and expression of the results | 11 |
| 9.1 Integration of peaks..... | 11 |
| 9.2 Calculation of inorganic arsenic in the samples | 11 |
| 10 Test report..... | 12 |
| Annex A (informative) Examples of typical chromatographic separation..... | 13 |
| Annex B (informative) Results of the inter-laboratory study..... | 17 |
| B.1 Inter-laboratory tests..... | 17 |
| B.2 Statistical results for determination of inorganic arsenic..... | 17 |
| Annex ZA (informative) Relationship of this European Standard and the essential requirements of Regulation (EU) 2019/1009 making available on the market of EU fertilising products aimed to be covered | 19 |
| Bibliography | 20 |

European foreword

This document (prEN 17706:2023) has been prepared by Technical Committee CEN/TC 455 “Plant Biostimulants”, the secretariat of which is held by AFNOR.

This document is currently submitted to the CEN enquiry.

This document will supersede CEN/TS 17706:2022.

This document has been prepared under a Standardization Request given to CEN by the European Commission and the European Free Trade Association, and supports essential requirements of EU Directive(s) / Regulation(s).

For relationship with EU Directive(s) / Regulation(s), see informative Annex ZA, which is an integral part of this document.

Introduction

This document was prepared by the experts of CEN/TC 455 “Plant Biostimulants”. The European Committee for Standardization (CEN) was requested by the European Commission (EC) to draft European standards or European standardization deliverables to support the implementation of Regulation (EU) 2019/1009 of 5 June 2019 laying down rules on the making available on the market of EU fertilizing products (“FPR” or “Fertilising Products Regulation”).

This standardization request, presented as M/564 and M/564 Amd1, also contributes to the Communication on “Innovating for Sustainable Growth: A Bio economy for Europe”. Working Group 4 “Other safety parameters” was created to develop a work programme as part of this request. Technical committee CEN/TC 455 “Plant Biostimulants” was established to carry out the work programme that will prepare a series of standards. The interest in biostimulants has increased significantly in Europe as a valuable tool to use in agriculture. Standardization was identified as having an important role in order to promote the use of biostimulants. The work of CEN/TC 455 seeks to improve the reliability of the supply chain, thereby boosting the confidence of farmers, industry, and consumers in biostimulants, and will promote and support commercialisation of the European biostimulant industry. The preparation of this document is based on a standardization request to CEN by the European Commission and the European Free Trade Association (Mandate M/564) concerning the modernization of methods of analysis of fertilizers in the framework of Regulation (EU) 2019/1009 of the European Parliament and of the Council.

This document describes a procedure of extraction and measurement for the determination of inorganic arsenic in plant biostimulants. The standard is based on a mild acid oxidative extraction of the arsenic species followed by liquid chromatography (HPLC or IC) coupled to the element-specific detector ICP-MS for the determination of the mass fraction of iAs.

WARNING — Persons using this document should be familiar with usual laboratory practice. This document does not purport to address all of the safety issues, if any, associated with its use. It is the responsibility of the user to establish appropriate health and safety practices and to ensure compliance with any national regulatory conditions.

IMPORTANT — It is absolutely essential that tests conducted according to this document are carried out by suitably trained staff.

1 Scope

This document specifies a method for extraction, separation, and determination of inorganic arsenic (iAs) in plant biostimulants using anion-exchange high performance liquid chromatography (HPLC) or ion chromatography (IC) coupled to ICP-MS. This document is applicable to the fertilizing product blends belonging to PFC 7 where the EU fertilising product plant biostimulant contained in the blend represents the highest % by mass in the blend. In case of equal shares, the user can apply either this or the standard(s) applicable to the other component product(s).

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 17704:—¹, *Plant biostimulants — Determination of dry matter*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

3.1

analyte

parameter to be determined

3.2

blank calibration solution

solution prepared in the same way as the calibration solution but leaving out the analytes

3.3

blank test solution

solution prepared in the same way as the test sample solution but omitting the test portion

3.4

calibration solution

solution used to calibrate the instrument, prepared from stock solutions by adding acids, buffer, reference element and salts as needed

3.5

stock solution

solution with accurately known analyte concentration(s), prepared from pure chemicals

4 Principle

This document describes a method for the determination of inorganic arsenic in plant biostimulants. Inorganic arsenic consists of arsenite As (III) and arsenate, As(V). A representative test portion of the sample is treated with a diluted nitric acid and hydrogen peroxide solution in a heated water bath. By this means the sample is solubilised, arsenic species are extracted into solution and As (III) is oxidized to As(V). The inorganic arsenic is selectively separated from other arsenic compounds using anion exchange HPLC (High Performance Liquid Chromatography) coupled online to the element-specific detector ICP-MS (Inductively Coupled Plasma Mass Spectrometer) for the determination of the mass fraction of the inorganic arsenic. External calibration with solvent matrix-matched standards is used for the

¹ Under preparation

quantification of the amount of the inorganic arsenic. Alternatively, IC (ion chromatography) coupled to ICP-MS can be used.

A preliminary determination of the total arsenic in aqua regia extracts by ICP-AES (EN 17701-1² and EN 17701-2³) could reduce the number of the samples where the determination of iAs is necessary because if the content of aqua regia (total) extractable arsenic is lower than the legislative limit for iAs then the determination of iAs is not necessary.

5 Reagents

When using a method of high sensitivity like ICP-MS and HPLC the control of the blank levels of water, acid and other reagents is very important. The reagents shall be of adequate purity and of recognized analytical grade. The concentration of arsenic species in the reagents and water used shall be negligible and low enough not to affect the results of the determination. Generally ultra-pure water from a purification system and nitric acid of minimum p.a. quality is recommended.

5.1 Water with an electrical conductivity not higher than 0,1 mS/m at 25 °C, having a resistivity greater than 18,2 MΩ·cm.

5.2 Nitric acid (HNO₃), concentrated, ≥ 65 % (mass fraction), mass concentration of approximately ρ (HNO₃) 1,4 g/ml.

Use only nitric acid available with high purity (minimum p.a. quality) in order to avoid potential contamination.

5.3 Hydrogen peroxide, H₂O₂ not less than 30 % (mass fraction).

High purity is essential to avoid potential contamination. Commercially available hydrogen peroxide for analysis should be tested for contamination of arsenic prior to use. It is necessary to prevent peroxide degradation and ensure the stability of the solution, this is in the discretion of the analyst to use only H₂O₂ of adequate quality

5.4 Extraction solution, 0,1 mol/l HNO₃ in 3 % (V/V) H₂O₂.

Add 6,5 ml of HNO₃ (5.2) and thereafter 100 ml of hydrogen peroxide (5.3) into 800 ml water (5.1) in a 1 000 ml volumetric flask. Fill the flask to the mark with water (5.1.). This solution is prepared on the day of use.

It is recommended that the total volume needed for the analysis is estimated and only this amount is produced in the day of use.

5.5 Ammonium carbonate, (NH₄)₂CO₃, mass fraction w ≥ 99,999 %, for preparation of the mobile phase solution.

5.6 Aqueous ammonia, (NH₃(aq).) w ≥ 25 %, for adjustment of pH in the mobile phase.

5.7 Methanol, (CH₃OH), HPLC grade, for preparation of the mobile phase solution.

5.8 Mobile phase, e.g. 50 mmol/l ammonium carbonate in 3 % methanol at pH 10,3.

Dissolve 4,80 g of ammonium carbonate (5.5) in approximately 800 ml water (5.1). Adjust the pH to 10,3 with aqueous ammonia (5.6) and add 30 ml of methanol (5.7) and then fill up to 1 000 ml with water (5.1). Prior to use filter the mobile phase solution through a 0,45 µm filter using a filtering device (6.4).

² Under preparation

³ Under preparation